

MBTタイプ



	Torq	Ang	R/L	INTERACTIVE .018		INTERACTIVE .022		PASSIVE .022	
				No Hook	With Hook	No Hook	With Hook	No Hook	With Hook
MAXILLA	1 1	+17	+4	R L	1475-18-1117 1475-18-2117		1485-22-1117 1485-22-2117		
	2 2	+10	+8	R L	1475-18-1210 1475-18-2210		1485-22-1210 1485-22-2210		
	3 3	0	+8	R L		1475-18-130HD 1475-18-230HD		1485-22-130HD 1485-22-230HD	1585-22-130HD 1585-22-230HD
	54 45	-7	0	R L	1475-18-1407 1475-18-2407	1475-18-1407HD 1475-18-2407HD	1485-22-1407 1485-22-2407	1485-22-1407HD 1485-22-2407HD	1585-22-1407 1585-22-2407
MANDIBULA	2 1 2	-6	0	R L	1475-18-4106 1475-18-3106		1485-22-4106 1485-22-3106		
	3 3	0	+3	R L		1475-18-430HD 1475-18-330HD		1485-22-430HD 1485-22-330HD	1585-22-430HD 1585-22-330HD
	4 4	-12	+2	R L	1475-18-4412 1475-18-3412	1475-18-4412HD 1475-18-3412HD	1485-22-4412 1485-22-3412	1485-22-4412HD 1485-22-3412HD	
	5 5	-17	+2	R L	1475-18-4517 1475-18-3517	1475-18-4517HD 1475-18-3517HD	1485-22-4517 1485-22-3517	1485-22-4517HD 1485-22-3517HD	

【症例セット】		INTERACTIVE .018		INTERACTIVE .022	
		No Hook	With Hook	No Hook	With Hook
UL5~5セット	MBTタイプ		1475MH		1485MH

*フックは3番のみ

ロスタイプ



	Torq	Ang	R/L	INTERACTIVE .018		INTERACTIVE .022	
				No Hook	With Hook	No Hook	With Hook
MAXILLA	1 1	+12	+5	R L	1475-18-1112 1475-18-2112		1485-22-1112 1485-22-2112
	2 2	+8	+9	R L	1475-18-128 1475-18-228		1485-22-128 1485-22-228
	3 3	-2	+9	R L		1475-18-1302HD 1475-18-2302HD	1485-22-1302HD 1485-22-2302HD
	54 45	-7	0	R L	1475-18-1407 1475-18-2407	1475-18-1407HD 1475-18-2407HD	1485-22-1407 1485-22-2407
MANDIBULA	2 1 2	-1	0	R L	1475-18-4101 1475-18-3101		1485-22-4101 1485-22-3101
	3 3	-11	+7	R L		1475-18-4311HD 1475-18-3311HD	1485-22-4311HD 1485-22-3311HD
	4 4	-17	0	R L	1475-18-4417 1475-18-3417	1475-18-4417HD 1475-18-3417HD	1485-22-4417 1485-22-3417
	5 5	-22	0	R L	1475-18-4522 1475-18-3522	1475-18-4522HD 1475-18-3522HD	1485-22-4522 1485-22-3522

ブラケットカラーコードはP.11参照

注:エンパワークリアのフックは遠心側についています

【症例セット】		INTERACTIVE .018		INTERACTIVE .022	
		No Hook	With Hook	No Hook	With Hook
UL5~5セット	ロスタイプ		1475RH		1485RH

*フックは3番のみ